

Instructions u-Slide

μ-Slide 8 Well high Bioinert



The ibidi product family is comprised of a variety of μ -Slides and μ -Dishes, which have all been designed for high-end microscopic analysis of fixed or living cells. The high optical quality of the material is similar to that of glass, so you can perform all kinds of fluorescence experiments with uncompromised resolution and choice of wavelength.

The μ -Slide 8 Well high Bioinert is an 8 well chamber slide in which cells can be cultivated and, subsequently, investigated with microscopical methods. This open μ -Slide (chambered coverslip) is intended for cell culture, immunofluorescence, live cell imaging, and high-end microscopy.

The Bioinert surface is a thin hydrogel layer that is covalently attached to the ibidi Polymer Coverslip No. 1.5. In contrast to standard ultra-low attachment (ULA) coatings, Bioinert is completely non-adherent and allows no binding of any biomolecule, even in long-term experiments. This makes Bioinert ideal for the culture and high resolution imaging of suspension cells and cell aggregates, like spheroids, organoids and embryoid bodies.

Material

ibidi μ -Slides, μ -Dishes, and μ -Plates are made of a polymer that has the highest optical quality. The polymer coverslip on the bottom exhibits extremely low birefringence and autofluorescence, similar to that of glass. Also, it is not possible to detach the bottom from the upper part. The μ -Slides, μ -Dishes, and μ -Plates are intended for one-time use and are not autoclavable, since they are only temperature-stable up to 80°C/175°F. Please note that gas exchange between the medium and the incubator's atmosphere occurs partially through the polymer coverslip, which should not be covered.

Optical Properties ibidi Polymer Coverslip

Refractive index n _D (589 nm)	1.52
Abbe number	56

Thickness No. 1.5 (180 μm)

Material Polymer coverslip

Please note! The ibidi Polymer Coverslip is compatible with certain types of immersion oil only. A list of suitable oils can be found on page 3.

Shipping and Storage

The μ -Slides, μ -Dishes and μ -Plates are sterilized and welded in a gas-permeable packaging. The shelf life under proper storage conditions (in a dry place, no direct sunlight) is listed in the following table.


	Conditions
Shipping conditions	Ambient
Storage conditions	RT (15-25°C)

	Shelf Life	
Bioinert	36 months	

Store the Bioinert products in a dry place (relative humidity <50%).

Geometry of the μ -Slide 8 Well high

The μ -Slide 8 Well ^{high} provides a standard slide format according to ISO 8037/1.



Geometry		
Outer dimensions in mm $(w \times l)$	25.5 × 75.5	
Number of wells	8	
Dimensions of wells in mm $(w \times l \times h)$	$9.4\times10.7\times9.3$	
Volume per well	300 µl	
Height with/without lid	10.8/9.5 mm	
Growth area per well	1.0 cm^2	
Coating area per well	$2.2\mathrm{cm}^2$	
Bottom	ibidi Polymer Coverslip	



u-Slide 8 Well high Bioinert

Instructions

Characteristics of the Bioinert Surface

Characteristics		
Bioinert surface thickness	200 nm	
Bioinert surface material	Polyol-based hydrogel	
Protein coatings	Not possible	

The Bioinert Surface

The Bioinert surface allows no adsorption, coating, or binding of proteins, antibodies, enzymes, and other biomolecules. Therefore, the Bioinert technology provides a stable passivation in cell-based assays for several days or even weeks. The hydrophilic Bioinert surface hinders any protein attachment, thus inhibiting subsequent cell attachment. The Bioinert surface is not biodegradable by cells allowing long-term assays with suspension cells and cell aggregates. Unlike with the ibiTreat and Uncoated surfaces, a coating is not possible.

Seeding Cells

Without a surface modification, Bioinert does not support direct cell adherence. Depending on your application the number of cells or cell aggregates might differ. Follow these steps for a general cell application protocol:

- Trypsinize and count cells as usual. Dilute the cell suspension to the desired concentration. Depending on your cell type, we recommend a 5-11 \times 10^4 cells/ml suspension.
- Apply 300 μl cell suspension into each well of the μ-Slide. Avoid shaking as this will result in inhomogeneous distribution of the cells.
- Cover with the supplied lid and incubate at 37°C and 5 % CO₂ as usual.

Undemanding cells can be left in their seeding medium for up to three days. However, best results might be achieved when the medium is exchanged every 1–2 days. To minimize cell disturbances, carefully aspirate $50\,\%$ of the old medium and replace it with fresh cell culture medium.

Tip:

Make sure to avoid uneven incubator shelves and microscope stages. Single cells or cell clusters might roll on one side over time. Please also avoid evaporation and temperature changes. Both will lead to convectional flow.

Microscopy

To analyze your cells, no special preparations are necessary. Cells can be directly observed live or fixed, preferably on an inverted microscope. The bottom cannot be removed. For optimal results in fluorescence microscopy and storage of fixed and stained samples, ibidi provides a mounting medium (50001) optimized for μ -Dishes, μ -Slides, and μ -Plates.

Chemical Compatibility

The following table provides some basic information on the chemical and solvent compatibility:

Chemical / Solvent	Compatibility
Methanol	yes
Ethanol	yes
Formaldehyde	yes
Acetone	yes, without lid
Mineral oil	no
Silicone oil	yes
Immersion oil	See Immersion Oil on page 3.



μ-Slide 8 Well high Bioinert

Instructions

Immersion Oil

When using oil immersion objectives with the ibidi Polymer Coverslip, use only the immersion oils specified in the table below. The use of any non-recommended oil could damage the ibidi Polymer Coverslip. The resulting leakage may harm objectives and microscope components. All immersion oils that are not listed in the table below should be considered as non-compatible.

Company	Product	Ordering No.	Lot Number	Test Date
ibidi	ibidi Immersion Oil	50101	16-12-27	01/2017
Cargille	Type A	16482	100592	01/2017
Cargille	Type HF	16245	92192	01/2017
Carl Roth	Immersion oil	X899.1	414220338	01/2017
Leica	Immersion Liquid	11513859	n.a.	03/2011
Nikon	Immersion Oil F2 30cc	MXA22192	n.a.	01/2020
Nikon	Silicone Immersion Oil 30cc	MXA22179	20191101	01/2020
Olympus	Silicone Immersion Oil	SIL300CS-30CC	N4190800	01/2017
Zeiss	Immersol 518 F	444960	160706	01/2017
Zeiss	Immersol W 2010	444969	101122	04/2012



Instructions

μ-Slide 8 Well high Bioinert

Ordering Information

u-Slide 8 Well high



Cat. No.	Description
80806	μ-Slide 8 Well high ibiTreat: #1.5 polymer coverslip, tissue culture treated, sterilized
80806-90	μ-Slide 8 Well ^{high} ibiTreat, Bulk Pack : #1.5 polymer coverslip, tissue culture treated, sterilized
80802	μ-Slide 8 Well high Collagen IV: #1.5 polymer coverslip, sterilized
80804	μ-Slide 8 Well high Poly-L-Lysine: #1.5 polymer coverslip, sterilized
80800	μ-Slide 8 Well high Bioinert : #1.5 polymer coverslip, surface passivation with Bioinert, sterilized
80801	μ-Slide 8 Well high Uncoated: #1.5 polymer coverslip, hydrophobic, sterilized
80807	$\mu\text{-Slide 8 Well}^{\ high}$ Glass Bottom: #1.5H (170 μm ± 5 $\mu m)$ D 263 M Schott glass, sterilized
80807-90	μ -Slide 8 Well ^{high} Glass Bottom, Bulk Pack: #1.5H (170 μm ±5 μm) D 263 M Schott glass, sterilized

μ-Slide 8 Well high Grid-500



Cat. No.	Description
80806-G500	μ -Slide 8 Well high ibiTreat Grid-500: #1.5 polymer coverslip, tissue culture treated, grid repeat distance 500 μm, sterilized

For research use only!

Further information can be found at ibidi.com. For questions and suggestions please contact us by e-mail *info@ibidi.de* or by telephone +49 (0)89/520 4617 0.

© ibidi GmbH, Lochhamer Schlag 11, 82166 Gräfelfing, Germany.