


ANIMA  LAB

iNSCOPIX

The Next Frontier of Discovery in Circuit Neuroscience



Proven Innovative Technology Trusted by Researchers Worldwide

Our complete neural circuit solutions
empower you to make ground-breaking discoveries
with speed and scientific rigor



Miniscope
Platforms
and
Applications



Data Analysis
Solutions



Enabling Success
with Expert Support

Cutting-edge technology

End-to-end integrated solutions

Over 150 publications in top-tier journals

Our Newest Miniscope System


nVue™ System: Dual color imaging during free behavior



Key Benefits

- Image two distinct neuronal populations simultaneously with single cell resolution
- Record two brain signals longitudinally over months
- Explore the intricacies of how two populations of neurons interact during free behavior

How it Works




nVue Applications

Poznaj wiele strategii obrazowania


Stacyjny + dynamiczny

Podwójnie dynamiczny



● wskaźnik statyczny lub dynamiczny ● wskaźnik dynamiczny
● Komórki ko-aktywne GCaMP i wskaźnik statyczny

Multicolor Registration Data



The Gold Standard Platforms for Imaging and Manipulating Brain Circuit Dynamics In Vivo

nVoke™ System



Integrated optogenetics and calcium imaging

Key Benefits

- Image large-scale brain activity with single-cell resolution during free behavior
- Record the same field of view over months with electronic focusing
- Integrate multimodal inputs with custom configurations and precise timing


nVista™ System



Calcium imaging during free behavior


nVoke Applications (simultaneous or sequential)

Obrazowanie jednej populacji komórek + opto-wzbudzenie/hamowanie projekcji z innej populacji komórek



○ Opsyna
● Wskaźnik wapnia

Obrazowanie + opto-wzbudzenie/hamowanie dwóch populacji komórek




Obrazowanie + opto-wzbudzenie/hamowanie tej samej populacji komórek




Image in Multiple Species

nVoke imaging and optogenetics in mouse




A, nVoke stimulation of BLA-NAC terminal projections plus simultaneous imaging of NAC neurons in a ChR2-cre;GCaMP6a;BLA-Cre mouse; **B**, BLA activation increases Ca^{2+} responses in NAC and is rewarding in mice exposed to place preference; **C**, ROIs identified (red, OG-stimulated; gray, non-stimulated) and **D**, data traces show stimulated (red) and non-stimulated (black) cells. Orange bars = mouse in OG-LED stimulation zone. Modified from [Stamatakis et al., \(2018\) Front Neurosci, 12:496](#).

nVista imaging in rat



A, Rat wearing nVista microscope with lens probe positioned in hippocampus CA1; **B**, Delta F/F image of cell activity in CA1; **C**, Spatial locations of individual cells color-coded and identified by PCA/ICA; **D**, Example Ca^{2+} activity traces ($n=6$; identified in color) from CA1 highlight individual neurons from panel C, over time (Inscopix data).

nVista imaging in marmoset



A, Common marmoset (*Callithrix jacchus*); **B**, Schematic of marmoset ladder climbing task; **C**, Image showing field of view (in primary motor cortex) through nVista (cells identified are color-coded in blue, red, and green); **D**, ROIs are color-coded to show cell firing sensitive to right (red), left (blue), or both (green) arm reach during climbing task. Images modified from [Kondo et al. \(2018\) Cell Reports 24, 2191-2195](#).


Workflow Solutions

Advanced workflow solutions to accelerate your imaging success

A Simplified Workflow

Ready-to-Image Viruses and ProView™ Integrated Lenses

ProView™ Express Probes




Ready-to-Image Viruses

- GCaMP viruses optimized for several brain regions in mice and rats
- Ensures consistent and reproducible GCaMP6 expression every time
- Available in ready-to-inject aliquots for immediate use

Mouse brain	Neocortex, hippocampus CA1, dorsal striatum
Virus construct	AAV1.Camk2a.GCaMP6f.WPRE.bGHpA
Rat brain	Prefrontal cortex, hippocampus CA1
Virus construct	AAV1.Camk2a.GCaMP6m.WPRE.bGHpA

ProView Integrated Lenses and Express Probes

- Gradient Refractive Index (GRIN) lens with attached baseplate enables efficient surgery and reproducible imaging across cohorts
- Integrated Lens reduces number of surgeries and Express Probes have a virus-coating reducing surgery to a *single step*
- Express Probe enables cell type-specific, reproducible imaging across cohorts and increases yield



Inscopix Commutator

- Reduces supervision during imaging sessions
- Responds to 360° of movement in any direction and reduces cable wear and tear



Software Solutions


Powerful software options to acquire and process Ca²⁺ imaging videos with Inscopix Data Acquisition (IDAS) and Data Processing Software (IDPS)

IDAS

Robust acquisition software makes Ca²⁺ imaging efficient

Key Benefits

- Live stream and record Ca²⁺ dynamics, crop and compress to speed up data transfer
- Digital focus enables multiplane imaging to maximize number of cells captured
- Advanced configurations allow flexible control of optogenetic stimulation and external devices



IDAS captures live Ca²⁺ dynamics in freely moving animals. *Multiplane imaging* enables the user to image up to three planes in a single session and maximize cell numbers imaged for more powerful results.

IDPS

Easily navigate Ca²⁺ imaging data using our complete streamlined processing pipeline; no scripting expertise is required

Minimize
data footprint



Motion
correction



Normalize
fluorescence



Cell
identification



Cell
registration




Event
detection



Key Benefits

- Automatically identify unique cells and extract Ca²⁺ traces using PCA/ICA or Inscopix CNMFe
- Perform cell registration: **longitudinal** – to identify the same cells over days, or **multicolor** – to identify two cell populations
- View and export high-resolution data into most formats for custom analyses and generate high-res manuscript images; simultaneously integrate Ca²⁺ videos with behavioral data



IDPS streamlines your raw Ca²⁺ videos with our intuitive processing pipeline to easily identify Ca²⁺ events. *Integrated visualization* enables viewing up to four panels of data (raw image, cell identification, live animal behavior and Ca²⁺ traces) for more impactful results.

Enabling Success with Expert Scientific Support

From experimental design to data analysis,
from grant to manuscript submission, we help
make great science happen faster

Field Scientific Consultants

Schedule virtual or
in person lab visits
with Inscopix's sales
representative Animalab



On-Demand Learning

Access online resources
including detailed
protocols, videos, and
user guides



Grant Assistance Program

Receive personalized
consultations on grant and
manuscript submissions for that
next neural circuit discovery



www.animalab.eu

info@animalab.eu

iNSCOPIX

Innovative Imaging to Empower Discovery in Brain Health

©2021 Inscopix Inc. The Inscopix logo, nVue, nVoke, nVista, and ProView are trademarks and/or registered trademarks of Inscopix.

Research use only. Not for use in diagnostic procedures.
MKC-MKT-004143 REV C

2462 Embarcadero Way
Palo Alto, CA 94303
888.874.6458
info@inscopix.com

www.inscopix.com