



INSTECH

F · S · T
 FINE SCIENCE TOOLS




Find the comprehensive offer of animal feeding needles for all rodent species which are the basic tool for feeding and intake drugs and different compounds types. Available models are made of polypropylene, stainless steel and polyurethane.

Plastic feeding tubes




- Flexible to reduce the risk of perforation
- Soft tip to minimize trauma
- Disposable – no cleaning and no risk of residual compound
- Faster to train new technicians
- Made of polypropylene (PP)



Stainless steel feeding tubes

- Reusable and autoclavable
- Easily bandable
- Available as curved and straight models
- Finished with a small ball to protect animal during feeding



Polyurethane feeding tubes

- Large inner diameter for easier delivery of viscous compounds
- Straight rounded tip
- Made of polyurethane (PU)
- One of the models offered is dedicated for dosing of gelatin capsules



Plastic feeding tubes

Barrel-tipped plastic feeding tubes				Packing / product codes				
Species	Lenght	Gauge	Ø OD	5 pcs.	50 pcs.	250 pcs.	Material	Brand
mouse pup	2.5 cm	22 ga	1.2 mm	-	FTP-22-25-50	FTP-22-25	PP	Instech
	3 cm	20 ga	1.6 mm	-	FTP-20-30-50	FTP-20-30	PP	Instech
mouse pup/adult	3.8 cm	20 ga	1.6 mm	-	FTP-20-38-50	FTP-20-38	PP	Instech
	3.8 cm	20 ga	1.8 mm	18961-20	-	-	PP	FST
mouse adult	3 cm	18 ga	2.3 mm	-	FTP-18-30-50	FTP-18-30	PP	Instech
	3.8 cm	16 ga	1.6 mm	-	FTPU-16-38-50	FTPU-16-38	PU	Instech
rat pup	3.8 cm	22 ga	1.2 mm	-	FTP-22-38-50	FTP-22-38	PP	Instech
rat 100-200 g	3.8 cm	18 ga	2.3 mm	-	FTP-18-38-50	FTP-18-38	PP	Instech
	3.8 cm	18 ga	2 mm	18961-18	-	-	PP	FST
	5 cm	18 ga	2.3 mm	-	FTP-18-50-50	FTP-18-50	PP	Instech
	6.4 cm	15 ga	3 mm	-	FTP-15-64-50	FTP-15-64	PP	Instech
rat ca. 200 g	5.2 cm	17 ga	2.4 mm	18961-17	-	-	PP	FST
	8.5 cm	17 ga	2.4 mm	18961-85	-	-	PP	FST
	7.5 cm	18 ga	2.3 mm	-	FTP-18-75-50	FTP-18-75	PP	Instech
	7.5 cm	18 ga	2 mm	18961-75	-	-	PP	FST
rat 200-300 g	8.5 cm	Capsule #9 doser	3.3 mm	-	FTPU-C9-85	-	PU	Instech
	8.8 cm	13 ga	2.4 mm	-	FTPU-13-88-50	FTPU-13-88	PU	Instech
rat 300+ g	7.8 cm	15 ga	3 mm	-	FTP-15-78-50	FTP-15-78	PP	Instech
	7.8 cm	15 ga	2.6 mm	18961-15	-	-	PP	FST
	10 cm	15 ga	3 mm	-	FTP-15-100-50	FTP-15-100	PP	Instech
large rodent	9 cm	13 ga	3.5 mm	-	FTP-13-90-50	FTP-13-90		Instech
	15 cm	13 ga	3.5 mm	-	FTP-13-150-50	FTP-13-150	PP	Instech
	38 cm	13 ga	-	-	FTPU-13-380-50	FTPU-13-380	PU	Instech

Stainless steel feeding tubes

Ball-tipped, stainless steel, feeding tubes				Product codes / packing			
Species	Lenght	Gauge	Ø OD	Straight models	Curved models	Packing	Brand
mouse pup	2.5 cm	24 ga	1.25 mm	18061-24	-	1 pcs.	FST
	2.5 cm	22 ga	1.2 mm	FTSS-22S-25		6 pcs.	Instech
	2.5 cm	22 ga	1.25 mm	18061-22	-	1 pcs.	FST
mouse adult	2.5 cm	20 ga	2.0 mm	QMI-01-25	-	6 pcs.	Qamar
	3.8 cm	20 ga	2.25 mm	18064-20	18065-20	1 pcs.	FST
	3.8 cm	20 ga	2.0 mm	QMI-02-38	-	6 pcs.	Qamar
rat 100-200 g	3.8 cm	18 ga	2.25 mm	18061-38	-	1 pcs.	FST
	5 cm	18 ga	2.25 mm	-	18061-50	1 pcs.	FST
rat 200-300 g	5.1 cm	18 ga	2.0 mm	QMI-03-51	-	6 pcs.	Qamar
	7.6 cm	18 ga	2.0 mm	QMI-04-76	-	6 pcs.	Qamar
rat 300+ g	7.5 cm	16 ga	3 mm	-	18061-75	1 pcs.	FST
	7.6 cm	16 ga	2.8 mm	QMI-05-76	-	6 pcs.	Qamar
	10 cm	16 ga	3 mm	-	18061-10	1 pcs.	FST