


**Your trusted global partner  
in life-saving research.**

For over 50 years, Allentown has provided critical solutions to the biomedical research community. Our high-quality equipment and advanced vivarium systems continue to define and exceed industry standards. Like everything we do, our unparalleled customer service is stamped with the Allentown hallmarks of integrity, dedication and care.


# Enrichment

FOR RODENTS



## OVERVIEW

Providing critical solutions with integrity and care.



## Richer Enrichment


The primary aim of environmental enrichment is to enhance animal well-being by providing animals with sensory and motor stimulation, through structures and resources that facilitate the expression of species-typical behaviors and promote psychological well-being through physical exercise, manipulative activities, and cognitive challenges according to species-specific characteristics. *Excerpted from The eighth edition of The Guide.*

**Elevated Enrichment** – Manufactured in sturdy polycarbonate, the Enrichment Loft provides the highest level of enrichment, privacy and safety for your rodent populations. Allentown offers both mouse and rat enrichment, compatible with the NexGen Mouse, NexGen MAX, NexGen Rat 900, NexGen Rat 1800 and 75 JAG cage systems.


**Cage Level Enrichment** – Manufactured in sturdy polycarbonate, Allentown offers the Mouse Mezzanine and the Safe Harbor Mouse Retreat. Both fit nicely into the cage without restricting access to water or food supplies while providing space, activity area and shelter. These devices work in a variety of Allentown and non-Allentown cage systems.

## OPTIONS


### Elevated Enrichment




**NexGen Mouse Enrichment Loft**  
60 mm x 78 mm x 102 mm | 45 g  
Item # 315680



**NexGen Rat 900 and Rat 1800 Enrichment Loft**  
189 mm x 106 mm x 90 mm | 136 g  
Item # 315584




**NexGen MAX Enrichment Loft**  
83 mm x 69 mm x 118 mm | 45 g  
Item # 316961



**75 JAG Enrichment Loft**  
145 mm x 64 mm x 63 mm | 45 g  
Item # 316415

### Cage Level Enrichment



**Safe Harbor Mouse Retreat™**  
64 mm x 70 mm x 114 mm | 50 g  
Item # SHMR001

Cage Level Enrichment works in a large variety of Allentown and non-Allentown cage types.



**Mouse Mezzanine™**  
160 mm x 70 mm x 170 mm | 86 g  
Item # SMMEZZ