


Part Number: **C10PU-MCA1459**



*not to scale*

**Full Description:** Catheter for mouse carotid artery, 1 to 2Fr, 8cm, square tip, collar @ 0.9cm. Fits 25ga

**Unit:** each      **Pouch Desc:** Mouse CAC PU 1-2Fr 8cm      **Pouch:** 3x7.5

**Approx Dead Vol:** 9  $\mu$ L      **Connects to:** 25ga      **EO Sterilized:**

**Intended Species:** mouse      **Vessel:** carotid artery

**Usage Notes:**

## TUBING

**Tip:** square      *Length 1*      *Length 2*  
**L1:** 2.5 cm      **L2:** 5.5 cm  
**Material (IDxOD):** PU 1Fr soft (.18x.35 $\pm$ .05mm)      **Material (IDxOD):** PU 2Fr (.43x.69 $\pm$ .03mm)

## FEATURES

<i>Distal (tip) end</i>		<i>Location</i>	<i>Fixed</i>
<b>Feature A:</b>	Small silicone collar for 1Fr (0.6mm OD), 1mm wide	a: 0.9 cm	<input type="checkbox"/>
<b>Feature B:</b>	Small silicone collar for 3Fr (1.3mm OD), 2mm wide	b: 2.5 cm	<input checked="" type="checkbox"/>
<b>Feature C:</b>		c: cm	<input type="checkbox"/>
<i>Proximal end</i>			
<b>Feature Z:</b>		z: cm	<input type="checkbox"/>
<b>Feature Y:</b>		y: cm	<input type="checkbox"/>
<b>Feature X:</b>		x: cm	<input type="checkbox"/>

## PARTS ADDED TO POUCH

**Part 1:**      **Part 2:**