



# RoVent<sup>®</sup> and RoVent<sup>®</sup> Jr.

Automated ventilator



RoVent and RoVent Jr. small animal automated ventilator are ideal for mice, rats and small animals up to 1,250 grams. Both are fully automated — simply enter the animal's weight and press run.

The Rovent includes the RightTemp temperature monitor and homeothermic warming module.

## Features & Benefits:

- **Fully automated with touchscreen control**  
*Safely ventilate using auto-calculated settings based on animal weight*
- **Volume control and pressure control**  
*Reduce footprint with one, dual control ventilator*
- **Multi-animal capability\***  
*Greater efficiency*

*\*Our clinical department will be pleased to review your study protocol to confirm that our product will meet your specific study requirements.*

**Standard 1-Year Parts & Labor Warranty** (Premium 3-year warranty available).  
Visit our website for more information.

Order your ventilator:

**888.572.8887**



**Kent Scientific**  
CORPORATION

Choose the ventilator that fits your research needs:

## Comparison Chart RoVent® & RoVent® Jr.



	RoVent Jr.	RoVent
Volume ventilation	●	●
Pressure ventilation	●	●
Pulse ventilation mode	○	●
PEEP	○	●
Assist mode	○	●
Sigh breath	○	●
I:E ratio	●	●
Tidal volume – measured	●	●
Auto set-up	Limited	●
User profiles	○	●
Alarms (parameter)	○	●
Trigger in	○	●
Analog out	○	●
Digital out	○	●
Sync out	○	●
ADI integration	○	●
Data output	○	●
Internal data logging	○	●
Realtime clock	○	●
Temperature monitor and homeothermic control module	○	Included (valued at \$1,000)

Order your ventilator:

**888.572.8887**  
kentscientific.com



**Kent Scientific**  
CORPORATION